

# *in vivo*-jetPEI®, an alternative to lipid-based reagents or viral vectors for nucleic acid-mediated therapies




Malik Hellal, Valérie Moreau, Fanny Prémartin, Thibaut Benchimol, Valérie Kedinger, Alengo Nyamay'antu, Pascale Belguise, Patrick Erbacher  
Polyplus-transfection®, Vectura, 75 rue Marguerite Perey, 67400 Illkirch, France

## Abstract


Nucleic acids have considerable potential as therapeutic agents in the treatment of pathologies including genetic diseases, viral infections, and cancer therapies. The major challenge for the use of nucleic acids in therapy lies in safe delivery of these anionic macromolecules to their intended sites of action. Our cationic polymer-based reagent, *in vivo*-jetPEI®, was tested for plasmid DNA and siRNA *in vivo* delivery through intravenous and intraperitoneal injections in mice to evaluate its efficacy, safety and biodistribution.

Here we demonstrate that *in vivo*-jetPEI® can efficiently and safely deliver various nucleic acids *in vivo* to target a wide range of tissues, through various routes of administrations.



Furthermore, *in vivo*-jetPEI® is widely acknowledged to deliver nucleic acids in animals; and coherently is selected as the delivery vector of choice in several drug development programs, notably for cancer therapies. To fulfill all the quality requirements associated to its use in Human, Polyplus® supplies a cGMP grade *in vivo*-jetPEI® reagent for a growing number of clinical trials.



## A wide variety of targeted organs depending on the administration route




## Safe method of delivery, with no major inflammatory response triggered




Bonnet et al., 2008

Bonnet, M.E., P. Erbacher, and A.L. Bolcato-Bellemin. 2008. *Pharmaceutical research*. 25:2972-2982.  
 Bonnet, M.E., J.B. Gossart, E. Benoit, M. Messmer, O. Zounib, V. Moreau, J.P. Behr, N. Lenne-Samuel, V. Kedinger, A. Meulle, P. Erbacher, and A.L. Bolcato-Bellemin. 2013. *Journal of controlled release*. 170:183-190.  
 Buscaill, L., B. Bournet, F. Vernejoul, G. Cambois, H. Lulka, N. Hanoun, M. Dufresne, A. Meulle, A. Vignolle-Vidoni, L. Ligat, N. Saint-Laurent, F. Pont, S. Dejean, M. Gayral, F. Martins, J. Torrisani, O. Barbey, F. Gross, R. Guimbaud, P. Otal, F. Lopez, G. Tiraby, and P. Cordelier. 2015. *Molecular therapy*. 23:779-789.  
 Matouk, I., E. Raveh, P. Ohana, R.A. Lail, E. Gershain, M. Gilon, N. De Groot, A. Czerniak, and A. Hochberg. 2013. *International journal of molecular sciences*. 14:4298-4316.  
 Sidi, A.A., P. Ohana, S. Benjamin, M. Shalev, J.H. Ransom, D. Lamm, A. Hochberg, and I. Leibovitch. 2008. *The Journal of urology*. 180:2379-2383.

## Dose response survival study of siRNA/*in vivo*-jetPEI® complexes




## *in vivo*-jetPEI® enables a seamless transition from Drug Discovery to Human clinical trials



*in vivo*-jetPEI® is available for applications ranging from *in vivo* Fundamental research or Proof-of-Concept experiments to Preclinical studies in animals (pharmacodynamics, biodistribution, toxicology studies...). Then, to meet quality and regulatory requirements for Human use in clinical trials and on the market, a cGMP-compliant grade is available: *in vivo*-jetPEI® GMP. In addition, Polyplus® also offers Regulatory support with:

- a Drug Master File (DMF) submitted to the US FDA that can be cross-referenced for IND applications & BLA,
- a Documentation describing the Chemistry, Manufacturing and Control (CMC) section for IMPD submission in Europe,
- a Quality agreement

## Clinical trials with *in vivo*-jetPEI® GMP



*in vivo*-jetPEI® has been selected as a nucleic acid delivery vector for the development of a growing number of nucleic acid-mediated therapies. Type of nucleic acid delivered, administration route and therapeutic application are very diverse.

## Conclusion

- ✦ **GMP compliant:** Manufactured in compliance with US and EU cGMP guidelines since 2007
- ✦ **Successful and proven:** trusted excipient used in several Human clinical trials worldwide (Buscaill et al., 2015; Sidi et al., 2008; Matouk et al., 2013) for different applications such as in cancer therapy, immunization, modulation of blood-brain barrier, etc.
- ✦ **Ready-to-use:** 2-step protocol requiring no equipment or formulation expertise
- ✦ **Polyvalent:** Used for *in vivo* delivery of any nucleic acid, to target any organ, in any animal model. Plus, different administration routes can be used, including systemic delivery
- ✦ **Fully supported:** Expert Regulatory and Scientific Support teams