

Technical specifications of TripleA 100 µl Tips

Applicable product: WB 5030



Eigenschaften:

Volume	100 µl
Low retention	Yes
Sterile	Yes
Filter	Yes
Packaging	Rack with hinge
Graduations	At 10, 50, and 100 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

The dimension and design of the 100 µl tip including the position of the filter. Dimensions are in millimeters.



SECTION A-A
SCALE: 3 : 1

