

Technical specifications TripleA 1250 µl Tips

Applicable products: WB 5074 - WB 5174 - WB 5570 - WB 5670

Characteristics:

Volume	1250 µl
Low retention	Yes
Sterile	Yes
Filter	Yes (WB 5074), Nee (WB 5174, WB 5570, WB 5670)
Packaging	Rack with hinge (WB 5074 , WB 5174), Bulk (WB 5570), Eco refill (WB 5670)
Graduations	At 100, 500, and 1000 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

The dimension and design of the 1250 µl tip including the position of the filter. Dimensions are in millimeters.

