

TripleA Tips - Pipettor Compatibility Chart

General Information

Based on in-house testing and information gathered from users. As such, it may be subject to change. In addition, the table lists only a limited selection of variable-volume, single-channel manual pipettes. To get 100% certainty that TripleA tips fit your pipettes, please contact us.

Brand	Model	Volume (µl)	TripleA catalog numbers end with:					
			10/15	20/30	40/45	50/55	64	70/74 /78
			10µl	20µl, 100µl	200µl	300µl	1000µl	1250µl
Biohit	mLine / Proline	0.1-3, 0.5-10	•					
	mLine / Proline	2-20, 5-50, 10-100, 20-200		•	•	•		
	mLine / Proline	100-1000					•	•
Brand	Transferpette	0.1-1, 0.1-2.5, 0.5-10, 2-20	•					
	Transferpette	2-20, 5-50, 10-100, 20-200		•	•	•		
	Transferpette	100-1000					•	•
Eppendorf	Research Plus / Reference	0.1-2.5, 0.5-10, 2-20 micro	•					
	Research Plus / Reference	2-20, 10-100, 20-200		•	•	•		
	Research Plus	30-300				•		
	Research Plus / Reference	100-1000					•	•
Thermo	Finnpette 4500	0.2-2, micro, 0.5-5, micro	•					
	Finnpette 4500	0.5-5, 2-20, 5-50, 10-100, 20-200		•	•	•		
	Finnpette 4500	100-1000					•	•
Gilson	Pipetman P2, P10 Classic/ G/L	0.2-2, 1-10	•					
	Pipetman P20, P100, P200 Classic/G/L	2-20, 20-100, 50-200		•	•	•		
	Pipetman P1000 Classic/G/L	200-1000					•	•
Oxford	Accupet Evo / Pro	0.1-2.5 , 0.5-10	•					
	Accupet Evo / Pro	2-20 , 5-50 , 10-100 , 20-200		•	•	•		
	Accupet Evo / Pro	100-1000					•	•
	Benchmate	0.2-2.0 , 0.5-10	•					
	Benchmate	2-20 , 10-100 , 20-200		•	•	•		
Rainin	Benchmate	100-1000					•	•
	Pipet-Lite SL- / Classic PR-	0.1-2, 0.5-10	•					
	Pipet-Lite SL- / Classic PR-	2-20, 10-100, 20-200		•	•	•		
	Pipet-Lite SL-	20-300				•		
Socorex	Pipet-Lite SL- / Classic PR-	100-1000					•	•
	Calibra 822 / Acura 825	0.2-2, 0.1-2, 1-10 micro, 0.5-10	•					
	Calibra 822 / Acura 825	1-10, 2-20, 5-50, 10-100, 20-200		•	•	•		
	Calibra 822 / Acura 825	100-1000					•	•