

Technical specifications of TripleA 300 µl Tips

Applicable products: WB 5050 - WB 5150 - WB 5155 - WB 5550 - WB 5650

Characteristics:

Volume	300 µl
Low retention	Yes
Sterile	Yes
Filter	Yes (WB 5040), No (WB 5150, WB 5155, WB 5550, WB 5650)
Packaging	Rack with hinge (WB 5050, WB 5150, WB 5155), Bulk (WB 5550), Eco refill (WB 5650)
Graduations	At 10, 50, 100 and 200 µl
Tip opening	Tapered for full sample release
Color	Colorless
Compatibility	i.a. Biohit, Brand, Eppendorf, Thermo/Finn, Gilson, Rainin, Socorex

For further information on product sterility and material composition please refer to the appendix

- COA Westburg tripleA tips
- Material composition sheet Westburg-plastics

The dimension and design of the 300 µl tip including the position of the filter. Dimensions are in millimeters.

