


DiluCell™




IMPLEN DiluCell™ ...

make it virtually dilute your sample!



DiluCell™ 10



- low sample volume necessary
- suitable for proteins and microbial cell cultures
- measurement possible starting at 340 nm
- automatic dilution of factor 10
- no risk of contamination
- significantly delayed cell sedimentation compared to standard cuvettes

| | |
|--------------------|---|
| Dimension | 12 mm x 12 mm x 45 mm |
| Pathlength | 1.0 mm (DC 10) (identifiable by special marking) |
| Automatic dilution | Factor 10 (DC 10) |
| Volume | DC 10: 300 µl max. sample volume |
| UV/VIS transparent | Starting from 340 nm |
| Availability | Approved for usage with the OD600 DiluPhotometer™ |

Features and specifications are subject to change without notice.



Implen GmbH
ISO 9001:2015



Europe, Asia, South Pacific, Middle East and Africa
Implen GmbH
Schatzbogen 52 | 81829 München | Germany
Tel. +49 89 7263718 20 | Fax +49 89 7263718 54 | info@implen.de | www.implen.de

North and South America
Implen, Inc.
31194 La Baya Drive, Unit 104 | Westlake Village | CA 91362 | USA
Tel. +1 818 748-6400 | Fax +1 818 449-0416 | info@implen.com | www.implen.com