

Anti-Digoxigenin (DIG) mab 11C8

Catalogue number	Package size
847-0102006701	100 µg
847-0102006703	1 mg

For research use only.

Product:	monoclonal antibody
Clone:	11C8
Host:	mouse
Isotype:	IgG1 _k
Immunogen:	Conjugate of bovine serum albumin and DIG
Specificity:	DIG
Purification:	purified by affinity chromatography
Purity:	> 95% by SDS-PAGE
Solution:	PBS pH 7.4 without additives
Storage:	at -20 °C (repeated thawing and freezing should be avoided)
Applications:	ELISA, Western Blot
Dilution:	1 µg/ml in Western Blot
Reactivity:	DIG, PCR products labelled with DIG, proteins labelled with DIG

